

# Experiment - 1

PAGE NO.: 001

DATE: 10 Aug

**Object:** To determine hardness of water by HCl.

**Apparatus:** Burette, Pipette, conical flask, beaker, volumetric flask, water sample, std N/40 HCl, methylorange.

## Procedure:

- \* Wash all apparatus with dist. water
- \* fill burette with N/40 HCl, note initial reading.
- \* Pipette out 10 ml of water sample into conical flask, add 2 drop of methylorange.
- \* Add HCl with constant string.
- \* Sharp colour change is noticed from Yellow to Cherry Red. Note the reading.
- \* Repeat till 2 concordant reading.
- \* Now, titrate with boil water in same manner.

## Observation Table:

S.No	Vol of water sample taken	Vol of N/40 HCl consumed				A = a - b = 2 - 0.4 = 1.6
		Before boiling	after boiling	concordant reading a	concordant reading b	
1	10 ml	2.0	0.4			
2	10 ml	2.0	0.4	2.0	0.4	
3	10 ml	2.0	0.4			

## Calculation:

Temporary hardness:  $1.6 \times \frac{1000}{10} \times \frac{50}{40} \text{ mg/l}$   
 $= 200 \text{ ppm}$

formula used =  $\frac{A \times 1000}{10 \text{ of N/40 HCl}}$

N  
24/8/24

## Experiment = 2

**Object:** To determine the hardness of water sample by complexometric method using EDTA

**Apparatus:** Burette, pipette, conical flask, beaker, burner, std  $M/100$  EDTA sol<sup>n</sup>, Eriochrome Black T indicator,  $NH_4Cl - NH_4OH$  buffer sol<sup>n</sup> of pH 10, hard sample.

**Procedure:** Pipette out 10 ml of std hard water in conical flask.

- \* Add 2 ml of buffer sol<sup>n</sup>, 2 drops of EBT as indicator.
- \* Titrate till colour changes wine red to Blue
- \* Repeat till concordant reading
- \* Pipette out 10 ml of boiled water, add 2 ml of buffer sol<sup>n</sup> and 2 drops of EBT
- \* Titrate with  $M/100$  EDTA sol<sup>n</sup> till blue colour appears.
- \* This corresponds to permanent hardness
- \* Same repeat with sample water.

**Observation Table:**

Titration = EDTA VS Std. hard water (10 ml)		Burette reading		Concordant Reading
S.No	Vol of water	Initial	Final	
1	10 ml	0	7	7
2	10 ml	0	7.5	
3	10 ml	0	7	

2. EDTA VIS Sample hard water ( $V_2$  ml)

S.No.	Vol of water	Burette reading		Concordant Reading
		Initial	Final	
1.	10 ml.	0	1.5	1.5
2.	10 ml	0	1.5	
3.	10 ml			

3. EDTA VIS Boiled water ( $V_3$  ml)

S.No.	Vol of boiled water	Burette reading		Concordant Reading
		Initial	Final	
1.	10 ml	0	0.5	0.5
2.	10 ml	0	0.5	
3.	10 ml.			

Calculation: Total hardness =  $1000 \times \frac{1.5}{7}$

Permanent hardness =  $1000 \times \frac{0.5}{7} \Rightarrow$

Temporary hardness =  $1000 \left( \frac{1.5 - 0.5}{7} \right) \Rightarrow 142.85 \text{ ppm}$

P  
20/8/24

Experiment = 3

**Object =** To determine the conductivity of given water samples by conductivity meter

**Apparatus and reagent required -**

Conductivity meter, water sample, distilled water, 0.01 N KCl solution, Beakers etc.

**Procedure =** \* Switch on the conductivity meter for about 15 to 20 min. for its warm up.

- \* Wash the conductivity cell in distilled water
- \* Set the function switch to check position
- \* Display must read 1.000 otherwise set it with CAL control knob of conductivity meter
- \* Adjust the temperature knob to room temperature
- \* Dip the washed and dried conductivity cell in the beaker containing sample solution
- \* After completion of experiment switch off the apparatus and put the cell in distilled water

**Observation**

$$l/a = 0.984$$

S. No	Water Sample Sol <sup>n</sup>	Conductivity
1	A	0.50 x 0.984
2	B	2.94 x 0.984
3	C	6.23 x 0.984

**Result =**

The conductivity of given water sample are as follow

Sample	Conductivity
A	0.492
B	2.89
C	6.130

$$K = \frac{l}{a} \times \text{conduct}$$

*[Signature]*  
31/8/24