

Experiment - 1

PAGE NO.: 01

DATE: 10 Aug

Object: To determine hardness of water by HCl.

Apparatus: Burette, Pipette, conical flask, beaker, volumetric flask, water sample, std N/40 HCl, methylorange

Procedure:

- ★ Wash all apparatus with dist. water
- ★ fill burette with N/40 HCl, note initial reading
- ★ Pipette out 10 ml of water sample into conical flask, add 2 drop of methylorange.
- ★ Add HCl with constant stirring.
- ★ Sharp colour change is noticed from Yellow to Cherry Red. Note the reading
- ★ Repeat till 2 concordant reading
- Now, titrate with boil water in same manner

Observation Table:

S.No	Vol of water sample taken	Vol of N/40 HCl consumed				A = a - b = 2 - 0.4 ⇒ 1.6
		Before boiling	after boiling	concordant reading a	concordant reading b	
1	10 ml	2.0	0.4			
2	10 ml	2.0	0.4	2.0	0.4	
3	10 ml	2.0	0.4			

Calculation:

Temporary hardness: $1.6 \times \frac{1000}{10} \times \frac{50}{40} \text{ mg/l}$

$= 200 \text{ ppm}$

formula used = $\frac{A \times 1000}{10 \text{ of N/40 HCl}}$

K

24/8/24

Experiment = 2

Object: To determine the hardness of water sample by complexometric method using EDTA.

Apparatus: Burette, pipette, conical flask, beaker, burner, std M/100 EDTA solⁿ, Eriochrome Black T indicator, $\text{NH}_4\text{Cl} - \text{NH}_4\text{OH}$ buffer solⁿ of pH 10, Hard sample.

Procedure: Pipette out 10 ml of std hard water in conical flask.

- * Add 1 ml of buffer solⁿ, 2 drops of EBT as indicator.
- * Titrate till colour changes wine red to Blue
- * Repeat till concordant reading.
- * Pipette out 10 ml of boiled water, add 2 ml of buffer solⁿ and 2 drops of EBT
- * Titrate with M/100 EDTA solⁿ till blue colour appears.
- * This corresponds to permanent hardness
- * Same repeat with sample water.

Observation Table:

1) Titration of EDTA VS Sld. Hard water (10 ml)		Burette reading		Concordant Reading
S.No	Vol of water	Initial	Final	
1.	10 ml	0	7	7
2.	10 ml	0	7.5	
3.	10 ml	0	7	

2. EDTA VIS Sample hard water (V_2 ml)

S.No.	Vol of water	Burette reading		Concordant Reading
		Initial	Final	
1.	10 ml.	0	1.5	1.5
2.	10 ml	0	1.5	
3.	10 ml			

3. EDTA VIS Boiled water (V_3 ml)

S.No.	Vol of boiled water	Burette reading		Concordant Reading
		Initial	Final	
1.	10 ml	0	0.5	0.5
2.	10 ml	0	0.5	
3.	10 ml.			

Calculation: Total hardness = $1000 \times \frac{1.5}{7}$

Permanent hardness = $1000 \times \frac{0.5}{7} \Rightarrow$

Temporary hardness = $1000 \left(\frac{1.5 - 0.5}{7} \right) \Rightarrow 142.85 \text{ ppm.}$